

## Alba PCR Grade Water

Cat. No.	Pack Size
MTE0054	10 x 2 ml
MTE0055	25 x 2 ml
MTE0056	5 x 10 ml

### Storage:

Store at 20-25°C for up to 12 months.

### Description

The **Alba PCR Grade H<sub>2</sub>O** is manufactured and tested for use as a component in PCR reactions.

This product is validated for use with highly sensitive PCR reaction including those with eubacterial primer sets. **Alba PCR Grade H<sub>2</sub>O** is stored in convenient 2 ml tubes to avoid unnecessary freeze-thaw of the whole water supply.

### Applications:

- Manufactured specifically for use in PCR reactions
- RNase, DNase and DNA contamination free
- Sterile and UV irradiated

### Quality Control:

Free of RNase, DNase and DNA contamination.

## Alba Magnesium Chloride Solution

Cat. No.	Pack Size
<b>MTE00578</b>	2 ml
<b>MTE0058</b>	5 x 2 ml

### Storage:

Store at -20°C for up to 12 months

### Description

**Alba MgCl<sub>2</sub> Solution** can be used for optimisation of magnesium ion concentration in PCR reactions.

### Quality Control:

Amplification of a single-copy gene from human genomic