

Vitrési Methods File Release Document

Methods File Version 4

Summary

The purpose of this document is to detail the updates incorporated into the Vitrési Methods File Version 4.

Features and Updates

Update	Description
Addition of new method program	The following kits can be performed on the SAS Vitrési when using Methods File version 4: SAS-1 Alk Phos Kit (REF:200801, 200802)
Updates to Pentavalent Screen (300406) methods	Updated 'Pentavalent Screen' and 'Urine and Serum Pentavalent Screen' methods to include back and forth steps between the EP bay and stainer to allow an extra blot and wash step. 'Pentavalent Screen' method name changed to 'Serum Pentavalent Screen' for applicator. 'Urine Pentavalent Screen' applicator method created to allow urine only method. Applicator methods numbered 1 and 3. Stainer methods numbered 2 and 4 with only a 'Pentavalent Screen' named method which covers serum, serum and urine and urine only methods.
Update to two-part kit method names - changed from reference of the top level code to the gel kit code.	The following have been changed: 200300 to 200301 300300 to 300301 102200 to 102201 104200 to 104201 201300 to 201301 301000 to 301001

Compatibility

Methods file version 4 is only compatible with Firmware version 7654.0100.01.00.

Upgrade

To upgrade the methods file, please refer to HL-2-P-3598 Vitrésí Methods File and Firmware Upgrade Instructions.

The Vitrésí Methods File can be found via the download link;

<https://www.helena-biosciences.com/vitresi-downloads>

Contact

For assistance and support, please contact Helena Biosciences Technical Support.